



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At5g27460 (At5g27460)

Product Code	CSB-BP664263DOA
Abbreviation	At5g27460
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E911
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLIRCVAVAK TRFNTTGVVS RLVSSLADGS DTSSVANRNS LKEILRKNGP RRSVTSLLEQ RIDSGHAVSL SELRLISKRL IRSNRYDLAL QMMEWMENQK DIEFSVYDIA LRLDLIKTH GLKQGEEYFE KLLHSSVSMR VAKSAYLPLL RAYVKNKMKVK EAEALMEKLN GLGFLVTPHP FNEMMKLYEA SGQYEKVVMM VSMMKGNKIP RNVLSYNLWM NACCEVSGVA AVETVYKEMV GDKSVEVGWS SLCTLANVYI KSGFDEKARL VLEDAEKMLN RSNRLGYFFL ITLYASLGNK EGVVRLWEVS KSVCGRISCV NYICVLSSLV KTGDLLEEAE VFSEWEAQCF NYDVRVSNVL LGAYVRNGEI RKAESLHGCV LERGGTPNYK TWEILMEGWV KCENMEKAID AMHQVFVLMR RCHWRPSHNI VMAIAEYFEK EEKIEEATAY VRDLHRLGLA SLPLYRLLLR MHEHAKRPAY DIYEMMKLKD L
Source	Baculovirus
Target Names	At5g27460
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At5g27460
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.