



Recombinant Arabidopsis thaliana Probable F-box protein At5g25300 (At5g25300)

Product Code	CSB-YP664264DOA
Abbreviation	At5g25300
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E960
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDTTNSSRWS ELPIDILRSL LEQKGGCAIY NPDEAKDYKT KRDSSGIRFL ANSVNWFLVL DSRSNLYIID VFSEKKIDLP PLESIKNGLY SLEQVGDEKK KEYIVVWFFD KGTEYIAFCK NGEDHYRDIP IRKNVCKELQ DLYDMLLHGG NNLYISKTRQ SIRKLRFFEE EGFIDVDNSE ILPFRKISFY IDGARFSNNI AVTTSREVLL VQNFFYETTR YRSFRLYKGD PNPDLNATIN NPYPMVELCS LGDVSASGFG DRKRSCLD MC MFNLVTKKLL RFPDLSNLKL TDARCTTTSS GNNWKFSIIV SEVIDSKFLE LRFADLDSES QQQTKLKFSM LKSEFVDMNK EDSQIEINEK ETKINQEHDQ SDETQAKRRR SLTSSVWQDL VSVGGEANVK EIDRSMETPN SSPPWSELPG DILRSVFERL SFVDFQRAKQ TCPIKRKSN CLRLWLITFG GLERLPHLSS VCFISWT
Source	Yeast
Target Names	At5g25300
Protein Names	Recommended name: Probable F-box protein At5g25300
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.