



# Recombinant Arabidopsis thaliana Probable aspartic protease At2g35615 (At2g35615)

<b>Product Code</b>	CSB-MP664270DOA
<b>Abbreviation</b>	At2g35615
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3EBM5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GHPKNFSVEL IHRDSPLSPI YNPQITVTDR LNA AFLRSVS RSRRFNHQLS QTDLQSG LIG ADGEFFMSIT IGTPPIKVFA IADTGSDLTW VQCKPCQQCY KENGPIFDKK KSSTYKSEPC DSRNCQALSS TERGCDESNN ICKYRYSYGD QSFSKGDVAT ETVSIDSASG SPVSFPGTVF GCGYNNGGTF DETGSGIIGL GGHLSLISQ LGSSISKKFS YCLSHKSATT NGTSVINLGT NSIPSSLSKD SGVVSTPLVD KEPLTYYYLT LEAISVGKKK IPYTGSSYNP NDDGILSETS GNIIDS GTT LTLLEAGFFD KFSSAVEESV TGAKRVSDPQ GLLSHCFKSG SAEIGLPEIT VHFTGADVRL SPINAFVKLS EDMVCLSMVP TTEVAIYGNF AQMDFLVGYD LETRTVSFQH MDCSANL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At2g35615
<b>Protein Names</b>	Recommended name: Probable aspartic protease At2g35615 EC= 3.4.-.-
<b>Expression Region</b>	21-447
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.