



Recombinant Arabidopsis thaliana F-box protein SKIP23 (SKIP23)

Product Code	CSB-MP664271DOA
Abbreviation	SKIP23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3EBZ2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVDWSTLPKD LLDLISKSL SSFDLIQFRS VCSSWRSAAE PKSPLPTHHL PILPDNGGSL FPDSAVGFRL SQRSILLIKP HEPCIESDSF GWLIKVEEDL NVPRKVTLLD PLCDTRNSIP ENFPRVLDMS KFKVRELGRE FKLHYFNTVG DIVESLYLEK AVVKYLD CDG DYKFVLLTIH VSGKLAVFRS WDRAWTVIND MPSPYDDVML FDGRFFAVDN NGRTVVVDYS SLKLTLVASP VFGGDKKFLI ESCGEMLLVD MYLSLEAVEG DPGFVEEIFE HPAFYMNERT VKFKVYRFVE REESWVDVYD LEDKMLFLGD DSTFSASASD ILPLCDGSSV FFNGNVFNGE DLGAMQDRDL GVFD FRSGKI ELVQKLPEYA KLFWPPPPWI TSHARAEDHI GETSSQS
Source	Mammalian cell
Target Names	SKIP23
Protein Names	Recommended name: F-box protein SKIP23 Alternative name(s): SKP1-interacting partner 23
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.