



Recombinant Arabidopsis thaliana Putative F-box protein At1g70960 (At1g70960)

Product Code	CSB-BP664272DOA
Abbreviation	At1g70960
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ECE2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MVNTSFETLP RHMQMEILSR VPLKFLMKFM CVSKKWASII RGEEFREDYL FQSMKRPRVL FVIDHREYLP IKPEAFFHSV YQEDQPLLLS GKQRMRTFET PLVQVFQPIR GLICQQGYGK IVICNPGLKK FRSLPQIKVH KGAPMRTFFG YDEDKDVFKV LCITWLRNGK RSEVSKEYLV YTMGSDEESS SWRLITCEHD HAPVTEGLFK GGVLYYGAKS NNGKSVMSF NVNSEDFSVI ELEVEISPYW RLVNYKGDIA LMNNIEDSLY HSREFEMWVR NEVTGNWDRT SIKISHWNGT VDGKTFYFKG TIGTKELAFA PDYWFGEQHF VSYYDTETKN LRRFDIEGMV DQDDFVRTFF DHVDSTWLI
Source	Baculovirus
Target Names	At1g70960
Protein Names	Recommended name: Putative F-box protein At1g70960
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.