



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At1g19930 (At1g19930)

Product Code	CSB-MP664275DOA
Abbreviation	At1g19930
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ED92
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKKNKVPTEL IFSLPNDLLV NILARVSRDL YPILSLVSKR FSSVLTLPPEL YQTRSLVGLT ENCLYVCLLS SRADRIPSWF KLCRRPILAS DTRKSSGYVL ATIPIPHSPP LHRSSLVAVG SNIYNIGGSI SQQSSSVSI LDCWSHTWLE GPSMQVEREY PSASLLDGKI YVTGGCRLTF HGCGDQTDNV VVDGKLHSCG GYKGVAYNPN DGRWDSLSE MNLGLKWSSS CVIENVIYY YHNENIKWYD TKVRSWRTL NGLKTLPRFAR YANVRLADYG GKMAFWDFK TGCGNQNRMI WCAIALERR NNEEIWGNVE WSDAVLPHQV PMAYVCEYVV AVNV
Source	Mammalian cell
Target Names	At1g19930
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At1g19930
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.