



# Recombinant Pseudoalteromonas haloplanktis Fructose-1,6-bisphosphatase class 1 1 (fbp1)

<b>Product Code</b>	CSB-BP664303PTP
<b>Abbreviation</b>	fbp1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3IFQ5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pseudoalteromonas haloplanktis (strain TAC 125)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRRLPPVLE DGCPRELISL IRTILAACKE ISFRVGQGEL SGVLGSTLDE NIQGETQKKL DVLTNQLLKD ILLESGYVKA IASEEEDYTV AGNPDAEYIV AFDPLDGSSN TDINSLVGTI FSVMKAPEGA DPADQSIFMQ PGINQVAAGY VLYGPSTILA LTTGKGTRFF TLDKTHGTFL LTQDFAKIPA DTNEYAINAS NQRHWQPAMQ NYINDLVAGD TGPRARNFNM RWIAAMVGDV HRVLSRGGLF TYPTDTKNPS QPNKLRLLYE ANPMAMLVEQ AGGIASGTGE RIMDIQPNAI HQRVAVILGS KNEVETCLGY HK
<b>Source</b>	Baculovirus
<b>Target Names</b>	fbp1
<b>Protein Names</b>	Recommended name: Fructose-1,6-bisphosphatase class 1 1 Short name= FBPase class 1 1 EC= 3.1.3.11 Alternative name(s): D-fructose-1,6- bisphosphate 1-phosphohydrolase class 1 1
<b>Expression Region</b>	1-322
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.