



Recombinant Bovine Recombining binding protein suppressor of hairless (RBPJ)

Product Code	CSB-BP664862BO
Abbreviation	RBPJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3SZ41
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MAPVVTGKFG ERPPKRLTR EAMRNYLKER GDQTVLILHA KVAQKSYGNE KRFFCPPPCV YLMGSGWKKK KEQMERDGCS EQESQPCAFI GIGNSDQEMQ QLNLEGKNYC TAKTLYISDS DKRKHFM LSV KMFYGN SDDI GVFLSKRIKV ISKPSKKKQS LKNADLCIAS GTKVALFNRL RSQTVSTRYL HVEGGNFHAS SQQWGAFYIH LLDDDESEGE EFTVRDGYIH YGQTVKLVCS VTGMALPRLI IRKVDKQTAL LDADDPVSQL HKCAFYLKDT ERMYLCLSQE RIIQFQATPC PKEPNKEMIN DGASWTIIST DKAEYTFYEG MGPVLAPVTP VPVVESLQLN GGGDVAMLEL TGQNFTPNLR VWFGDVEAET MYRCGESMLC VVPDISAFRE GWRWVRQP VQ VPVTLVRNDG IYSTSLTFT YTPEPGPRPH CSAAGAILRA NSSQVPPNES NTNSEGSYTN VSTNSTSVTS STATVVS
Source	Baculovirus
Target Names	RBPJ
Protein Names	Recommended name: Recombining binding protein suppressor of hairless Alternative name(s): J kappa-recombination signal-binding protein RBP-J kappa
Expression Region	1-487
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.