



Recombinant Mouse Protein FAM98A (Fam98a)

Product Code	CSB-EP664908MO-B
Abbreviation	Fam98a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3TJZ6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MECDLMETDI LESLEDLGYK GPLLDDGALL QAVSAGAASP EFTKLCAWLV SELRVLCKLE ENVQATNSPS EAEEFQLEVS GLLGEMNCPY PSLTSGDVTK RLLVEKNCLL LLTYLISELE AARMLCVNAP PKKAQEGGGS EVFQELKGIC IALGMSKPPA NITMFQFFSG IEKCLKETLA KVPPNHVGKP LLKKPMGPAH WEKIEAINQA IANEYEVRRK LLIKRLDVTV QSGWSDRAK SQTDKLAKVY QPKRSLSPK GKVSVAHLLA ARQDLSKILR TSSGSIREKT ACAINKVLMG RVPDRGGRPN EIEPPPEMP PWQKRQDGPQ QQAGGRGGGR GGYEHSSYGG RGGHEQGGRG GRGSYDHGGR GGGRGKNHQG GWTGGSGSG GGYQDGAYRD SGFQPGGYHG GHSGGYQAGG YGGFQTSSYT GSGYQGGGYQ QDNRYQDGGH HGERGSGRGG RGGRGGRGGR GSQGGGWGGR GSQTYHQGGQ FEQHFQHGYY QYSHSGFGQG RHYTS
Source	E.coli
Target Names	Fam98a
Protein Names	Recommended name: Protein FAM98A
Expression Region	1-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.