



Recombinant Mouse Coiled-coil domain-containing protein 102A (Ccdc102a)

Product Code	CSB-MP664910MO
Abbreviation	Ccdc102a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3TMW1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSHGSPRLA ESPQLSKSSL LTILGSPSPE RMGPADSLPP TPPSGTPSPG PPPSLPLTAP ALLADGDWES REELRLRELE EARARALQME KTMRWWSRCT ANWREKWSKV RAERNRAREE VRQLRQLDT LTKELAGARR ERQEAQGECE ARGRELARLR GVHSAADKTH DGPEPEREQE PVRDIGAERP PGSQELDLVE SLLKSRPEEP EGCWDACSVA AGGPRVSSGR QDRNRLPWED TASTEEDASK LTALRLRLDE SQKVLLKERE DKLALSKNIE KLEGELSQWK IKYEELSKTK QEMLKQLSIL KETHQDELGR MSEDLEDELG ARSSMDRKMA ELRGEMERLQ AENAAEWGRR ERLETEKLGL ERENKKLRAQ VGDLEEALAR RRRQNASALD CDLRASQAAL FEKNKELADL KHVHGKLLKQ FQEKVAELAH ANRRVEQHET EVKKLRLRVE ELKKELAQAE DELDEAHNQA RKLQRSLEDEQ TEQSENLQVQ LEHLQSRLRR QQQNAPLFGK IRSTRFGTEE AGDGASDLDE DEDLQIQVA
Source	Mammalian cell
Target Names	Ccdc102a
Protein Names	Recommended name: Coiled-coil domain-containing protein 102A
Expression Region	1-549
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.