



# Recombinant Mouse E3 ubiquitin-protein ligase DTX3L (Dtx3l)

<b>Product Code</b>	CSB-EP664947MO
<b>Abbreviation</b>	Dtx3l
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3UIR3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ASSPDPSP LLVRLRESIP KAHRKLEIYF QSRASGGGEC SVQPVGPSAP DTYEVKFLKK ADKEKVLKKS EHEMLVHNKP VTIVLETTKK PVEDLRPRLP SLTQPVETPS SRPPSLTGSL DEALCDDIHP QDGLVSNSVD SVVQKIFLAV TAEINCDLLS KEQRASITTV CPHIIKSMEG SDGIKKVCGN FKDIEKIHFF LSEQLLEREQ KRKGSEQKRK CAPQKHTPPD VEREPDQSS IQVPVLLLEY FKHVNPGRLE FIEYKFGVNI EIQASSPNMV TVGFTSSPFG NVEEASQSFV RDFQKCSQSL KQDCISLEEH QRAKEVRQEL SRCFPKLLIK GQGRTLTLG SPADISAATE KVSQGLGLRP VKITASGYTT GIEVDSTRFK LLEPELLQEI SEIEQKFNTR GKVQEKGQKT CILFVPKDKD LDLSVQSYTG FTDAFQRATW QLRTEVLSLK GLGKERARLH NTKFADNFKK EHPNVHFVTS QESVTLTGLP HHLAQAMQYV SKRMGLAPSS GEKLAMDQET PMEISSSDPH GDQQENAALP APRGTSSSPA ASKGTEDYCV ICMDTISNKH VLPKCKHEFC TSCISKAMLI KPVCPVCLTS YGIQKGNQPE GTMSYSTQKG SLPGYEGCGT IVINYEIKDG IQTKEHPNPG KAYHGTRRTA YLPDNTEGRK VLDLLHEAFK HRLTFTIGYS RATGVSDVIT WNDIHHKTSK FGGPANFGYP DPDYLRVKE ELKAKGIE
<b>Source</b>	E.coli
<b>Target Names</b>	Dtx3l
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase DTX3L EC= 6.3.2.- Alternative name(s): Protein deltex-3-like
<b>Expression Region</b>	2-748
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.