



Recombinant *Saccharomyces cerevisiae* RNA polymerase II transcription factor B subunit 5 (TFB5)

Product Code	CSB-BP665574SVG
Abbreviation	TFB5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E7C1
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MARARKGALV QCDPSIKALI LQIDAKMSDI VLEELDDTHL LVNPSKVEFV KHELNRLLSK NIYNPMDEEE NQ
Source	Baculovirus
Target Names	TFB5
Protein Names	Recommended name: RNA polymerase II transcription factor B subunit 5 Alternative name(s): General transcription and DNA repair factor IIH subunit TFB5 Short name= TFIIH subunit TFB5
Expression Region	1-72
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.