



# Recombinant Arabidopsis thaliana B3 domain-containing protein REM2 (REM2)

<b>Product Code</b>	CSB-MP665589DOA
<b>Abbreviation</b>	REM2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3E9T2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDDPAISSPT NKHFFVIVLH GHKSSPMIPA EFFSTYVEGK NHQSTKLKLT SDAFDRTWEV KLNRRRFAGG WENFSTVHSL QDDDVVIFRE IGDMTFHVTA SGRSFCKLQY ISSDEGRADE DDAWNVYSDS KNIVSKKKPR TEAESSKKNY YLVAHVTPSN LSRDRLYLIK AFARSNGLNE RCCEIDLMNK HGKSWTLGLR YSTVNDQASI SGGWRDFCLV NGLEAGSFYR FKLQVNGERP LLRLCYDTIP EEDCYKGNRK VNLSESASMK QGKFLTMTFK PYMLKSGQLR LSRPFAIENE IKEAEDITLV DKNNGNWPSY VASGDGQGGK YYLAKGWRFS CAANRLKTGE TFTLEFVRGE GTTPMLKFCS KSKDKIEQVP FDEGGETQEE TETRVKKRAR VSAERGSSRR NQASNKPNAD PQNLQPGQPH CSISDQVQKV KQIIVDTLTG IKRIRSELER KGQKLEASIL EIDALGEKVS EINKILK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	REM2
<b>Protein Names</b>	Recommended name: B3 domain-containing protein REM2 Alternative name(s): Protein REPRODUCTIVE MERISTEM 2
<b>Expression Region</b>	1-487
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.