



Recombinant Human Keratin-associated protein 13-3 (KRTAP13-3)

Product Code	CSB-BP666167HU
Abbreviation	KRTAP13-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3SY46
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSYNCCSRNF SSCSHGGYLH YPGSSCGSSY PSNLVYSTDL CSPSTCQLGS SLYRGCQETC WRPNSCQTLC VESSPCHTSC YYPRTHMLCN SCLTMHVGSR GFGSNSCCSL SCGSRSCSSL GCGSNGFRYL NYRIHTSPSQ SYRSRFCHPI YFPPRRWFHS SCYQPFCRSG FY
Source	Baculovirus
Target Names	KRTAP13-3
Protein Names	Recommended name: Keratin-associated protein 13-3
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.