



Recombinant Mouse ETS-related transcription factor Elf-3 (Elf3)

Product Code	CSB-YP666275MO
Abbreviation	Elf3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UPW2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAATCEISNV FSNYFNAMYS SEDPTLAPAP PTTFGTEDLV LTLNNQQMTL EGPGPQTRSQ RDRTDPLAVL HLAEKASWTS ERPQFWSKTQ VLEWISYQVE KNKYDASSID FSRCDMDGAT LCSCALEELR LVFGPLGDQL HAQLRDLTSN SSDELSWIIE LLEKDGMSFQ ESLGDSGPFQ QGSPFAQELL DDGRQASPHY CSTYGGGAPS PGSSDVSTAR TATPQSSHAS DSGGSDVDLD LTESKVFPRD GFPDYKKGEP KHGKRKRGRP RKLSKEYWDC LEGKSKHAP RGTHLWEFIR DILIHPELNE GLMKWENRHE GVFKFLRSEA VAQLWGQKKK NSNMTYEKLS RAMRYYYKRE ILERVDGRRL VYKFGKNSSG WKEEEVGESR N
Source	Yeast
Target Names	Elf3
Protein Names	Recommended name: ETS-related transcription factor Elf-3 Alternative name(s): E74-like factor 3 Epithelial-restricted with serine box Epithelium-restricted Ets protein ESX Epithelium-specific Ets transcription factor 1 Short name=
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.