



Recombinant Human Polycomb group RING finger protein 3 (PCGF3)

Product Code	CSB-BP667288HU
Abbreviation	PCGF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KNV8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLTRKIKLWD INAHITCRLC SGYLIDATTV TECLHTFCRS CLVKYLEENN TCPTCRIVIH QSHPLQYIGH DRTMQDIVYK LVPGLQEAEM RKQREFYHKL GMEVPGDIKG ETCSAQHLD SHRNGETKAD DSSNKEAAEE KPEEDNDYHR SDEQVSICLE CNSSKLRGLK RKWIRCSAQA TVLHLKKFIA KKLNLSSFNE LDILCNEEIL GKDHTLKFVV VTRWRFFKAP LLLHYRPKMD LL
Source	Baculovirus
Target Names	PCGF3
Protein Names	Recommended name: Polycomb group RING finger protein 3 Alternative name(s): RING finger protein 3A
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein contains a C3HC4 type RING finger, which is a motif known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.