



Recombinant Human Probable tubulin polyglutamylase TTLL9 (TTLL9)

Product Code	CSB-MP667484HU
Abbreviation	TTLL9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3SXZ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVPSREALLG PGTTAIRCPK KLQNNQNYKGGH GLSKGKEREQ RASIRFKTTL MNTLMDVLRH RPGWVEVKDE GEWDFYWCDV SWLRENFHDHT YMDEHVRISH FRNHVELTRK NYMVKNLKRK RKQLEREAGK LEAAKCDFFP KTFEMPCEYH LFVEEFRKNP GITWIMKPVA RSQGGIFLFL RRLKDIVDWR KDTRSSDDQK DDIPVENYVA QRYIENPYLI GGRKFDLRVY VLVMSVFAEC LLWSGHRQD VHLTNVAVQK TSPDYHPKKG CKWTLQRFRQ YLASKHGPEA VETLFRDIDN IFVKSLQSVQ KVIISDKHCF ELYGYDILID QDLKPWLLEV NASPSLTASS QEDYELKTCL LEDTLHVVDV EARLTGREKR VGGFDLMWND GPVSREEGAP DLSGMGNFVT NTHLGCVNDR KKQLRQLFCS LQVQKKASS
Source	Mammalian cell
Target Names	TTLL9
Protein Names	Recommended name: Probable tubulin polyglutamylase TTLL9 EC= 6.-.-. Alternative name(s): Tubulin--tyrosine ligase-like protein 9
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.