



Recombinant Mouse Coiled-coil domain-containing protein 61 (Ccdc61)

Product Code	CSB-MP667577MO
Abbreviation	Ccdc61
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UJV1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MDQPAGLQVD YIFRGVEHAV RVVVSGQVLE LEVEDRMTAD QWRGEFDANF IEDLTHKTGN FKQFSIFCNM LESALTQSSE SVTLDLLTYT DLESLRSRKL GGRPGALAPR SAQLNSKRYL ILIYSVEFDR IHYPLPLPYQ GKPDPVVLQG IIRSLKEELG HLRGLNGGQD ARETEIWHLR EQVTRLTSEK RELEAQLGRS REEALAGRAA RQEAESLRGL VRGLELELRQ ERGLGGRAAG RRSQDSRRLA KELEEVKASE RNLRARLKTL NSELAMYRRG RRTLPPVAVR EGRASSSRER STSRGRAVTR SSSRESSRGA RGRGRPAHPS PSPTGSRAPR FDPTAFVKAK EKKQREIRMK QQQQQRNRMG SGGSGDGPVS SWSQQTRPAA AVTGRGDVN RSRNRSSSVD SFRSRCSSVS SCSELEDFSH SVSRSRRCRG RGKPPSPTPW SRSKTKSTTQ ERSDHQRHLA SSGGWVPIKE YSSDFQGADM AEIDARLKAL QEYMNRLDTR S
Source	Mammalian cell
Target Names	Ccdc61
Protein Names	Recommended name: Coiled-coil domain-containing protein 61
Expression Region	1-511
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.