



Recombinant *Gossypium hirsutum* Polygalacturonase (G9)

Product Code	CSB-MP667733GHB
Abbreviation	G9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q39786
Product Type	Recombinant Protein
Immunogen Species	<i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>)
Purity	≥85% (SDS-PAGE)
Sequence	DAFD VVAKFGAKAD GKTDLSKPFL DAWKEACASV TPSTVVIPKG TYLLSKVNLE GPCKAPIEIN VQGTIQAPAD PSAFKDPNWV RFYSVENFKM FGGGIFDGQG SIAYEKNTCE NREFRSKLPV NIRFDLFTNA LIQDITSKDS KLFHINVFAC KNITLERLKI EAPDESPNTD GIHMGKSEGV NIIASDIKTG DDCISIGDGT KNMVIKEITC GPGHGISIGS LGKFQNEEPV EGIKISNCTI TNTSNGARIK TWPGEHGGAV SEIHFEDITM NNVSSPILID QQYCPWNKCK KNEESKVKLS NISFKNIRGT SALPEAIKFI CSGSSPCQNV ELADIDIKHN GAEPATSQCL NVKPITSGKL NPIPESGPVP KTPSATA
Source	Mammalian cell
Target Names	G9
Protein Names	Recommended name: Polygalacturonase Short name= PG EC= 3.2.1.15 Alternative name(s): Pectinase
Expression Region	27-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.