



Recombinant *Saccharomyces cerevisiae* rRNA-processing protein EFG1 (EFG1)

Product Code	CSB-EP668170SVG-B
Abbreviation	EFG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E705
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAKLQRKRSK ALGSSLEMSQ IMDAGTNKIK RRIRDLERLL KKKKDILPST VIIKERNLQ ALRLELQNN E LKNKIKANAK KYHMVRF FEK KKALRKYNRL LKKIKESGAD DKDLQQKLRA TKIELCYVIN FPKTEKYIAL YPNDTPSTDP KGVELTNLRR EQFLKLVAER MDANTLNVSF EEILKGKKLD EDSIGLTLSP DKDHEDGSQV SPTQDRKELD QVVGEDKDD FFE
Source	E.coli
Target Names	EFG1
Protein Names	Recommended name: rRNA-processing protein EFG1 Alternative name(s): Exit from G1 protein 1
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.