



Recombinant *Saccharomyces cerevisiae* Cytochrome c oxidase subunit 6B-like protein YMR244C-A (YMR244C-A)

Product Code	CSB-MP668188SVG
Abbreviation	YMR244C-A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3E846
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGLFSFDGGK KESQPPNTRS QRKLCWESRD AFFQCLDKAD ILDAMDPKNS KSIKSHCKVE NEKFEENCAH SWIKYFKEKR VIDFKREQTI KRIEQEAKQR ERNQ
Source	Mammalian cell
Target Names	COA6
Protein Names	Recommended name: Cytochrome c oxidase subunit 6B-like protein YMR244C-A
Expression Region	1-104
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.