



# Recombinant Human Putative acyl-coenzyme A thioesterase 6 (ACOT6)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP668210HU  |
| <b>Abbreviation</b>      | ACOT6   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q3I5F7  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MLQHPKVKGP SIALLGFSKG GDLCLSMASF LKGITATVLI NACVANTVAP<br>LHYKDMIIPK LVDDLKVKI TKSGFLTFMD TWSNPLEEHN HQSLVPLEKA<br>QVPFLFIVGM DDQSWKSEFY AQIASERLQA HGKERPQIIC YPETGHCIDP<br>PYFPPSRASV HAVLGEAIFY GGEPKAHKA QVDAWQQIQT FFHKHLNGKK<br>SVKHSKI   |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | ACOT6   |
| <b>Protein Names</b>     | Recommended name: Putative acyl-coenzyme A thioesterase 6 Short name=<br>Acyl-CoA thioesterase 6 EC= 3.1.2.-  |
| <b>Expression Region</b> | 1-207   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |