



Recombinant Rat Kelch domain-containing protein 2 (Klhdc2)

Product Code	CSB-MP668596RA
Abbreviation	Klhdc2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3KRE6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MADGNEDARA EDMPGPAFEG YEAMELACPA ERSGHVAVSD GRHMFVWGGY KSNQVRGLYD FYLPREELWI YNMETGRWKK INTEGDVPPS MSGSCAVCVD RVLVYFGGHH SRGNTNKFYM LDSRSADRVL QWERIDCQGI PPSSKDKLGV WVYKNKLIFF GGYGYLPEDK VLGTFEFD SFWNSSHPRG WNDHVHILDT EFAWSQPIT TGKAPSPRAA HACATVGNKG FVFGGRYRDA RMNDLHYLNL DTWEWNELIP QGICPVGRSW HSLTPVSSDH LFLFGGFTTE KQPLSDAWTY CISKNEWIQF NHPYVEKPRL WHTACASDEG EVIVFGGCAN NLLVHHRAAH SNEVLIFSVQ PKSLVRLSLE AVICFKEMLA NSWSCLPKHL LHSVNRFRGS NNTSGS
Source	Mammalian cell
Target Names	Klhdc2
Protein Names	Recommended name: Kelch domain-containing protein 2
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.