



Recombinant Rat ADP-ribosylation factor GTPase-activating protein 2 (Arfgap2)

Product Code	CSB-EP668675RA-B
Abbreviation	Arfgap2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3MID3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	AAGPSKSEI QTLFKRLRAI PTNKACFDCG AKSPSWASIT YGVFLCIDCS GVHRSLGVHL SFIRSTELDS NWSWLQLRCM QVGGNANATA FFRQHGCLAN DANTKYNSRA AQMYREKIRQ LGSTALARHG TDLWIDNMNS APSHSPEKGD SDFTEHTQA PAWDTAATDP SGTQQPALPS ESSSLAQPEP GPNTDLLGTS PQASLELKSS IIGKKKPAIAA KKGLGAKKGL GAQKVSNSQSF TEIERQAQVA EKLREQQAAD AKKQAEESMV ASMRLAYQEL QIDRKKEEEK LQNLEGKKRE QAERLGMGLV SRSSISHSVL SEMQMIEQET PLSAKSSRSQ LDFDDEVGTF ASGPPKYKDN PFSLGETFGS RWDSDAAWGM DRVEEKEPEV TISSIRPISE RTTSRREVES RISGLESSEA RQKFAGAKAI SSDMFFGREV DSEYEARSRL QQLSGSSAIS SSDLFGDVDG AHGGGTVSLG NVLPTADIAQ FKQGKSVAG KMAVLANGVM NSLQDRYGSY
Source	E.coli
Target Names	Arfgap2
Protein Names	Recommended name: ADP-ribosylation factor GTPase-activating protein 2 Short name= ARF GAP 2 Alternative name(s): GTPase-activating protein ZNF289 Zinc finger protein 289
Expression Region	2-520
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.