



# Recombinant Mouse Cullin-4A (Cul4a)

<b>Product Code</b>	CSB-EP668835MO
<b>Abbreviation</b>	Cul4a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3TCH7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADEGPRKGS VSALMGR TNG LTKPAALAGG PAKPGGTGGS RKLVIKNFRD RPRLPDNYTQ DTWRKLHEAV KAIQSSTSIR YNLEELYQAV ENLCSHKVSP TLYKQLRQVC EDHVQAQILP FREDSLDSVL FLKKINTCWQ DHCRQMIMIR SIFLFLDRTY VLQNSMLPSI WDMGLELFRN HIISDRMVQS KTIDGILLI GRERSGEAVD RLLRSLLSM LSDLQVYKDS FELKFLEETN CLYAAEGQRL MQDREVPEYL NHVSKRLEEE ADRVITYLDH STQKPLIACV EKQLLGEHLT AILQKGLEHL LDENRVPDLT QMYQLFSRVK GGQHALLQHW SEYIKTFGTT IVINPEKDKD MVQDLLDFKD KVDHVVEVCF QRNERFINLM KESFETFINK RPNKPAELIA KHVDSKLRAG NKEATDEELE RILDKIMILF RFIHGKDVFE AFYKKDLAKR LLVGKSASVD AEKSMLSKLK HECGAAFTSK LEGMFKDMEL SKDIMVHFKQ HMQNSAPGP IDLTVNILTM GYWPTYTPME VHLPPPEMVR QEVFKTFYLG KHSGRKLQWQ TTLGHAVLKA DFKEGKKEFQ VSLFQTLVLL MFNEGDGFSF EEIKMATGIE DSELRRTLQS LACGKARVLI KSPKGKEVED GDKFIFNADF KHKLFRIKIN QIQMKETVEE QVSTTERVFQ DRQYQIDAAI VRIMKMRKTL GHNLLVSELY NQLKFPVKPG DLKKRIESLI DRDYMERDKD SPNQYHYVA
<b>Source</b>	E.coli
<b>Target Names</b>	Cul4a
<b>Protein Names</b>	Recommended name: Cullin-4A Short name= CUL-4A
<b>Expression Region</b>	1-759
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .