



Recombinant Mouse Cullin-4A (Cul4a)

Product Code	CSB-MP668835MO
Abbreviation	Cul4a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3TCH7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MADEGPRKGS VSALMGR TNG LTKPAALAGG PAKPGGTGGS RKLVIKNFRD RPRLPDNYTQ DTWRKLHEAV KAIQSSTSIR YNLEELYQAV ENLCSHKVSP TLYKQLRQVC EDHVQAQILP FREDSLDSVL FLKKINTCWQ DHCRQMIMIR SIFLFLDRTY VLQNSMLPSI WDMGLELFRN HIISDRMVQS KTIDGILLI GRERSGEAVD RLLRSLLSM LSDLQVYKDS FELKFLEETN CLYAAEGQRL MQDREVPEYL NHVSKRLEEE ADRVITYLDH STQKPLIACV EKQLLGEHLT AILQKGLEHL LDENRVPDLT QMYQLFSRVK GGQHALLQHW SEYIKTFGTT IVINPEKDKD MVQDLLDFKD KVDHVVEVCF QRNERFINLM KESFETFINK RPNKPAELIA KHVDSKLRAG NKEATDEELE RILDKIMILF RFIHGKDVFE AFYKKDLAKR LLVGKSASVD AEKSMLSKLK HECGAAFTSK LEGMFKDMEL SKDIMVHFKQ HMQNSAPGP IDLTVNILTM GYWPTYTPME VHLPPPEMVR QEVFKTFYLG KHSGRKLQWQ TTLGHAVLKA DFKEGKKEFQ VSLFQTLVLL MFNEGDGFSF EEIKMATGIE DSELRRTLQS LACGKARVLI KSPKGKEVED GDKFIFNADF KHKLFRIKIN QIQMKETVEE QVSTTERVFQ DRQYQIDAAI VRIMKMRKTL GHNLLVSELY NQLKFPVKPG DLKKRIESLI DRDYMERDKD SPNQYHYVA
Source	Mammalian cell
Target Names	Cul4a
Protein Names	Recommended name: Cullin-4A Short name= CUL-4A
Expression Region	1-759
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.