



Recombinant Mouse Keratin, type I cytoskeletal 26 (Krt26)

Product Code	CSB-EP668844MO-B
Abbreviation	Krt26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3TRJ4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSFRFSGRSR VCSRAGSVRL SRGGAGFVAG NVCVGPGAEG SFSCTLGGLS SGGSFASEGS GRGSSIGFLN NEPGLFSGNE KVAMQNLNDR LALYLNHVSA LEEANTDLEK KIEGWYEKCG PGRGRRLDHD CSRYFSVIED LKRQILSMITT CNANLVLQND NARLTADDFK MKYENELALH QSVEADTNGL RRVLDELTLT TTDLEIQREA LSEELTYLQK NHEEEMVVLQ NASGGNINVE MNAAPSVDLT AMLNNMRAEY EDLAEQNRKD AEASFKEKSA SLQQQISDDA GAATAARNEL MELKRNLQTL EIELQSITAM KQSYENSLAE TEGNYAQLQ QIQEQIGARE EQLQQIRTET EGQKLEHEQL LGIKTCLEKE IDTYCNLLDG EEQRSESTSY KPKDGKPASE FNDSAETFA RTVAEELDQL GNLLSLRVHS VEEKSSKISN ITMEQRVPSK AP
Source	E.coli
Target Names	Krt26
Protein Names	Recommended name: Keratin, type I cytoskeletal 26 Alternative name(s): Cytokeratin-26 Short name= CK-26 Keratin-26 Short name= K26 Type I inner root sheath-specific keratin-K25irs2
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the keratin superfamily. This keratin protein is a type I keratin that is specific for the inner root sheath of the hair follicle. This gene exists in a cluster with other keratin genes on chromosome 17q21.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.