



Recombinant Mouse Protein zyg-11 homolog B (Zyg11b)

Product Code	CSB-BP668868MO
Abbreviation	Zyg11b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UFS0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MPEDQAHAAM EEASPYSLLD ICLSFLTTLNLEKFC SARQDGLTCLQEPGVF PQEVADRLLQ TIAFHGLLND GTVGIFRGNQ MRLKRACIRK AKISAVAFRK AFCHHKLVEL DATGVNADIT ITDIISGLGS NKWIIQQNLQC LVLNSLTLSL EDPYERCFSR LSGLRALSIT NVLFYNEDLA EVASLPRLES LDISNTSITD ITALLACKDR LKSLTMHHLK CLKMTTQIL DVVRELKHLN HLDISDDKQF TSDIALRLLQ QKDILPNLVS LDVSGRKHVT DKAVEAFIQQ RPSMQFVGLL ATDAGYSEFL MGKGHKLVSG EANETQIAEA LRRYSERAFV VREALFHLFS LTHVMEKTKP DILKLVVTGM RNHPMNLVQ LAASACVFNL TKQDLALGMP VRLADVTHL LLKAMEHFPN HQQLQKNCLL SLCSDRILQD VPFNRFEEAK LVMQWLCNHE DQNMQRMAVA IISILAAKLS TEQTAQLGAE LFIVRQLLQI VKQKTNQNSV DTTLKFTLSA LWNLDESPT TCRHFIEHQG LELFMRVLES FPTESIQQK VLGLLNNIAE VQELHSELMW KDFIDHISSL LHSVEVEVSY FAAGIIAHLI SRGEQAWTLS RSQRNSLLDD LHSAILKWPT PECEMVAYRS FNPFFPLLGC FTTPGVQLWA VWAMQHVCSK NPSRYCSMLI EEGGLQHLYN IKEHEQTDQPY VQQIAVAILD SLEKHIVRHG RPPPCKKQPQ ARLN
Source	Baculovirus
Target Names	Zyg11b
Protein Names	Recommended name: Protein zyg-11 homolog B
Expression Region	1-744
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.