



Recombinant Mouse Protein zyg-11 homolog B (Zyg11b)

Product Code	CSB-YP668868MO
Abbreviation	Zyg11b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UFS0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MPEDQAHAAM EEASPYSLLD ICLSFLTTLNLEKFC SARQDGLTCLQEPGVFPQEVADRLLQTIAFHGLLND GTVGIFRGNQMRLKRACIRK AKISAVAFRKAFC HHKLVEL DATGVNADIT ITDIISGLGS NKWIIQQNLQCLVLSLTLSEDPYERCFSLSGLRALSIT NVLFYNEDLA EVASLPRLS LDISNTSITDITALLACKDR LKSLTMHHLK CLKMTTQIL DVVRELKHLN HLDISDDKQFTSDIALRLLQKDILPNLVS LDVSGRKHVT DKAVEAFIQQ RPSMQFVGLLATDAGYSEFL MGKGHLKVSG EANETQIAEA LRRYSERAFF VREALFHLFSLTHVMEKTKP DILKLVVTGM RNHPMNLVQ LAASACVFNL TKQDLALGMPVRLADVTHL LKAMEHFVN HQQLQKNCLL SLCSDRILQD VPFNRFEEAKLVMQWLCNHE DQNMQRMAVA IISILAAKLS TEQTAQLGAE LFIVRQLLQIVKQKTNQNSV DTTLKFTLSA LWNLDESPT TCRHFIEHQG LELFMRVLESFPTESSIQQK VLGLLNIAE VQELHSELMW KDFIDHISSL LHSVEVEVSYFAAGIIAHLI SRGEQAWTLS RSQRNSLLDD LHSAILKWPT PECEMVAYRSFNPFFPLLGC FTTPGVQLWA VWAMQHVCSK NPSRYCSMLI EEGGLQHLYN IKEHEQTDQPY VQQIAVAILD SLEKHIVRHG RPPPCKKQPQ ARLN
Source	Yeast
Target Names	Zyg11b
Protein Names	Recommended name: Protein zyg-11 homolog B
Expression Region	1-744
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.