



Recombinant Mouse Polycomb group RING finger protein 5 (Pcgf5)

Product Code	CSB-YP668878MO
Abbreviation	Pcgf5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UK78
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MATQRKHLVK DFNPYITCYI CKGYLIKPTT VTECLHTFCK TCIVQHFEDS NDCPRCGNQV HETNPLEMLR LDNTLEEIIF KLVPGFREQE LQRELEFWKK NKPQENGQDD ISKVDKSKAD EEGDENQDDK DYHRSDPQIA ICLDCLRNGG QSGDNVVKGL MKKFIRCSTR VTVGTIKKFL SLKCLKPSSY ELDVLCNGEI MGKDHTMEFI YMTRWRLRGE NLCCCLNCSAS QVCSQDGTLY QSYPMVLQYR PRIDFG
Source	Yeast
Target Names	Pcgf5
Protein Names	Recommended name: Polycomb group RING finger protein 5 Alternative name(s): RING finger protein 159
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.