



Recombinant Bovine Protein BEX3 (NGFRAP1)

Product Code	CSB-MP669136BO
Abbreviation	NGFRAP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ZBJ6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MANIQENEE MEQPVQNGEE DRPLGGGEGH QPERNHRRGQ ARRLAPNFRW VIPNRQVNDG MGGEGDDMEI FMEEMREIRR KLRELQLRNC LRILMGELSN HHDHHDEFCL MP
Source	Mammalian cell
Target Names	BEX3
Protein Names	Recommended name: Protein BEX3 Alternative name(s): Brain-expressed X-linked protein 3 homolog Nerve growth factor receptor-associated protein 1 p75NTR-associated cell death executor
Expression Region	1-112
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.