



Recombinant Arabidopsis thaliana Oxygen-evolving enhancer protein 2-1, chloroplastic (PSBP1)

Product Code	CSB-MP669235DOA
Abbreviation	PSBP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q42029
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AYG EAANVFGKPK TNTDFLPYNG DGFKVQVPAK WNPSKEIEYP GQVLR FEDNF DATSNLNMV TPTDKKSITD YGSPEEFLSQ VNYLLGKQAY FGETASEGGF DNNAVATANI LESSSQEVGG KPYYYSVLT RTADGDEGGK HQLITATVNG GKLYICKAQA GDKRWFKGAR KFVESAATSF SVA
Source	Mammalian cell
Target Names	PSBP1
Protein Names	Recommended name: Oxygen-evolving enhancer protein 2-1, chloroplastic Short name= OEE2 Alternative name(s): 23 kDa subunit of oxygen evolving system of photosystem II 23 kDa thylakoid membrane protein OEC 23 kDa subunit Short na
Expression Region	78-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.