



Recombinant Arabidopsis thaliana Isopentenyl-diphosphate Delta-isomerase II, chloroplastic (IPP2)

Product Code	CSB-BP669244DOA
Abbreviation	IPP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q42553
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	AFSGT AMTDTKDAGM DAVQRRLMFE DECILVDETD RVVGHDSKYN CHLMENIEAK NLLHRAFSVF LFNSKYELL QRSNTKVTF PLVWTNTCCS HPLYRESELI QDNALGVRNA AQRKLLDELG IVAEDVPVDE FTPLGRMLYK APSDGKWGEH ELDYLLFIVR DVKVQPNPDE VAEIKYVSRE ELKELVKKAD AGEEGLKLSP WFRLVVDNFL MKWWDHVEKG TLVEAIDMKT IHKL
Source	Baculovirus
Target Names	IPP2
Protein Names	Recommended name: Isopentenyl-diphosphate Delta-isomerase II, chloroplastic EC= 5.3.3.2 Alternative name(s): Isopentenyl pyrophosphate isomerase II Short name= IPP isomerase II
Expression Region	46-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.