



Recombinant Citrullus lanatus Cysteine synthase

Product Code	CSB-EP669274CKU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q43317
Product Type	Recombinant Protein
Immunogen Species	Citrullus lanatus (Watermelon) (Citrullus vulgaris)
Purity	≥85% (SDS-PAGE)
Sequence	MADAKSTIAK DVTELIGHTP LVYLN RVVDG CVARVAAKLE MMEPCSSVKD RIGYSMISDA ENKGLITPGE SVLIEPTSGN TGIGLAFIAA AKGYRLIICM PASMSLERRT ILRAFGAELV LTDPARGMKG AVQKAEEIKA KTPNSYILQQ FENPANPKIH YETTGPEIWR GSGGKIDALV SGIGTGGTVT GAGKYLKEQN PNIKLYGVEP VESAILSGGK PGPHKIQQIG AGFIPGVLDV NLLDEVIQVS SEESIETAKL LALKEGLLVG ISSGAAAAAA IRIAKRPENA GKLIVAVFPS FGERYLSTVL FESVKRETEN MVFEP
Source	E.coli
Protein Names	Recommended name: Cysteine synthase Short name= CSase EC= 2.5.1.47 Alternative name(s): Beta-PA/CSase EC= 2.5.1.51 Beta-pyrazolylalanine synthase L-mimosine synthase EC= 2.5.1.52 O-acetylserine (thiol)-lyase
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.