



Recombinant Outer membrane lipoprotein blc

Product Code	CSB-YP669354DTO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46036
Product Type	Recombinant Protein
Immunogen Species	Citrobacter freundii
Purity	>85% (SDS-PAGE)
Sequence	CS SPTPPKGVTV VNNFDAKRYL GTWYEIARFD HRFERGLDKV TATYSLRDDG GINVINKGYN PDREMWQKTE GKAYFTGDPS TAALKVSFFG PFYGGYNVIA LDREYRHALV CGPDRDYLWI LSRTPTISDE MKQQMLAIAT REGFEVNKLI WVKQPGA
Source	Yeast
Protein Names	Recommended name: Outer membrane lipoprotein blc
Expression Region	19-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.