



Recombinant Beta-lactamase Toho-1 (bla)

| | |
|--------------------------|---|
| Product Code | CSB-EP669494ENL |
| Abbreviation | bla |
| Storage | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |
| Uniprot No. | Q47066 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli |
| Purity | >85% (SDS-PAGE) |
| Sequence | A NSVQQQLEAL EKSSGGRLGV ALINTADNSQ ILYRADERFA MCSTSKVMAA AAVLKQSESD KHLNQRVEI KKS DLVNYNP IAEKHVNGTM TLAELGAAAL QYSDNTAMNK LIAHLGGPDK VTAFARSLGD ETFRLDRTEP TLNTAIPGDP RDTTTP LAMA QTLKNLTLGK ALAETQRAQL VTWLKGNTTG SASIRAGLPK SWVVGDKTGS GDYGTNDIA VIW PENHAPL VLVTYFTQPE QKAERRRDIL AAAAKIVTHG F |
| Source | E.coli |
| Target Names | bla |
| Protein Names | Recommended name: Beta-lactamase Toho-1 EC= 3.5.2.6 |
| Expression Region | 30-291 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |