



Recombinant Polysialic acid biosynthesis protein P7 (neuC)

Product Code	CSB-BP669512ENL
Abbreviation	neuC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q47400
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	KKILYVTGS RAEYGIVRRL LTMLRETPEI QLDLAVTGMH CDNAYGNTIH IIEQDNFNII KVV DININTT SHTHILHSMS VCLNSFGDFF SNNTYDAVMV LGDRYEIFSV AIAASMHNIP LIHIGGEKT LANYDEFIRH SITKMSKLHL TSTEEYKKRV IQLGEKPGSV FNIGSLGAEN ALSLHLPNKQ ELELKYGSLL KRYFVVVFHP ETLSTQSVND QIDELLSAIS FFKNTHDFIF IGSNADTGSD IIQRKVKYFC KEYKFRYLIS IRSEDYLAMI KYSCGLIGNS SSGLIEVPSL KVATINIGDR QKGRVARGASV IDVPVEKNAI VRGINISQDE KFISVVQSSS NPYFKENALI NAVRIIKDFI KSKNKDYKDF YDIPECTTSY D
Source	Baculovirus
Target Names	neuC
Protein Names	Recommended name: Polysialic acid biosynthesis protein P7
Expression Region	2-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.