



# Recombinant Mouse Polycomb complex protein BMI-1 (Bmi1)

<b>Product Code</b>	CSB-EP002726MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P25916
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MHRTTRIKIT ELNPHLMCVL CGGYFIDATT IIECLHSFCK TCIVRYLETS KYCPICDVQV HKTRPLLIR SDKTLQDIVY KLVPGFLFKNE MKRRRDFYAA HPSADAANGS NEDRGEVADE EKRIITDDEI ISLSIEFFDQ SRLDRKVNKE KPKEEVNDKR YLRCPAAMTV MHLRKFLRSK MDIPNTFQID VMYEEEEPLKD YYTLMDIAYI YTWRRNGPLP LKYRVRPTCK RMKMSHQ RDG LTNAGELESD SGSDKANSPA GGVPSTSSCL PSPSTPVQSP HPQFP HISST MNGTSPNSPSA NHQSSFASRP RKSSLNGSSA TSSG
<b>Source</b>	E.coli
<b>Target Names</b>	Bmi1
<b>Protein Names</b>	Recommended name: Polycomb complex protein BMI-1 Alternative name(s): Polycomb group RING finger protein 4
<b>Expression Region</b>	1-324
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.