



Recombinant Human Coiled-coil domain-containing protein 37 (CCDC37)

Product Code	CSB-BP669847HU
Abbreviation	CCDC37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q494V2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSEIPSTIVS KNMTNDKNSL ESMNIISSSSS TEENPKKQAR KNEEHGPDPS ANPFHLSGDV DFFLLRDQER NKALSERQQQ KTMRVHQKMT YSSKVS AKHT SLRRQLQLED KQEDLEARAE AEHQRAFRDY TTWKLTLTKE KNVEPENMSG YIKQKRQMFL LQYALDVKRR EIQRLETLAT KEEARLERAE KSLEKDAALF DEFVRENDCS SVQAMRAAEK ETKAKIEKIL EIRDLTTQIV NIKSEISRFE DTLKHYKVYK DFLYKLSPE WLEEQEKHS FLKKAKEVSE ASKESSVNST PGDKGPGIKG KASSMWAKEG QGTKKPWRFL QTMRLGRSPS YLSSPQQGSQ PSESSGGDSR GSNSPIPTQ EDTSDGEEP QLYFTEPQQL LDVFRELEEQ NLSLIQNSQE TEKLEELSH TLKHTQIRMD REVNQLKQWV TTMMMSITKE EDTAAEELK ARVFHFGEYK GDQQDKLLES LNCKVLDVYR HCTGTQQEAN LGTVQMLTII EHQLDELLEN LEHVPQVKIE QAERAKEKER RIRLREEKLQ MQKILQEEHL QRARARAQAE IKKKRGRTL V CRSRPPAHRI KQQSEHTLMD KEEEEELLFFF T
Source	Baculovirus
Target Names	CFAP100
Protein Names	Recommended name: Coiled-coil domain-containing protein 37
Expression Region	1-611
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.