



# Recombinant Human Coiled-coil domain-containing protein 37 (CCDC37)

<b>Product Code</b>	CSB-MP669847HU
<b>Abbreviation</b>	CCDC37
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q494V2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSEIPSTIVS KNMTNDKNSL ESMNISSSSS TEENPKKQAR KNEEHGPDPS ANPFHLSGDV DFFLLRDQER NKALSERQQQ KTMRVHQKMT YSSKVS AKHT SLRRQLQLED KQEDLEARAE AEHQRAF RDY TTWKLTLTKE KNVEPENMSG YIKQKRQMFL LQYALDVKRR EIQRLETLAT KEEARLERAE KSLEKDAALF DEFVRENDCS SVQAMRAAEK ETKAKIEKIL EIRDLTTQIV NIKSEISRFE DTLKHYKVYK DFLYKLSPE WLEEQEKKHS FLKKAKEVSE ASKESSVNST PGDKGPGIKG KASSMWAKEG QGTKKPWRFL QTMRLGRSPS YLSSPQQGSQ PSESSGGDSR GSNSPIPTQ EDTSDGEEP QLYFTEPQQL LDVFRELEEQ NLSLIQNSQE TEKTLEELSH TLKHTQIRMD REVNQLKQWV TTMMSITKE EDTAAEELK ARVFHFGEYK GDQQDKLLES LNCKVLDVYR HCTGTQQEAN LGTVQMLTII EHQLDELLEN LEHVPQVKIE QAERAKEKER RIRLREEKLQ MQKILQEEHL QRARARAQAE IKKKRGRTL V CRSRPPAHRI KQQSEHTLMD KEEEEELLFFF T
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CFAP100
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 37
<b>Expression Region</b>	1-611
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.