



Recombinant Human Ankyrin repeat and death domain-containing protein 1A (ANKDD1A)

Product Code	CSB-MP669849HU
Abbreviation	ANKDD1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q495B1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQEELAWETD GLLPLERQLH EAARQNNVGR MQELIGRRVN TRARNHVGRV ALHWAAGAGH EQAVRLLLEH EAAVDEEDAV GALTEARLCF GMNALLLSAW FGHLRILQIL VNSGAKIHCE SKDGLTLLHC AAQKGHPVL AFIMEDLEDV ALDHVDKLGR TAFHRAAEHG QLDALDFLVG SGCDHNVKDK EGNTALHLAA GRGHMAVLQR LVDIGLDLEE QNAEGLTALH SAAGGSHPCD VQLLLRAGST VNALTQKNLS CLHYAALSGS EDVSRVLIHA GGCANVVDHQ GASPLHLAVR HNFPAVRLR INSDSDVNAV DNRQQTPLHL AAEHAWQDIA DMLLIAGVDL NLRDKQGKTA LAVAVRSNHV SLVDMIIKAD RFYRWEKDHP SDPSGKSLSF KQDHRQETQQ LRSVLWRLAS RYLQPREWKK LAYSWEFTEA HVDAIEQQWT GTRSYQEHGH RMLLIWLHGV ATAGENPSKA LFEGLVLAIGR RDLAGWSTMA RSQLTATSAS RVQMILVPQP PE
Source	Mammalian cell
Target Names	ANKDD1A
Protein Names	Recommended name: Ankyrin repeat and death domain-containing protein 1A
Expression Region	1-522
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.