



Recombinant Human mTERF domain-containing protein 3, mitochondrial (MTERFD3)

Product Code	CSB-EP669874HU-B
Abbreviation	MTERFD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q49AM1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	TTDKQ SSKENTRTVE KLYKCSVDIR KIRRLKGWVL LEDETYVEEI ANILQELGAD ETAVASILER CPEAIVCSPT AVNTQRKLWQ LVCKNEEELI KLIEQFPESF FTIKDQENQK LNVQFFQELG LKNVVISRLL TAAPNVFHNP VEKNKQMVRI LQESYLDVGG SEANMKVWLL KLLSQNPFIL LNSPTAIKET LEFLQEQGFT SFEILQLLSK LKGFLFQLCP RSIQNSISFS KNAFKCTDHD LKQLVLKCPA LLYYSVPVLE ERMQGLLREG ISIAQIRETP MVLELTPQIV QYRIRKLNSS GYRIKDGHLA NLNGSKKEFE ANFGKIQAKK VRPLFNPVAP LNVEE
Source	E.coli
Target Names	MTERF2
Protein Names	Recommended name: mTERF domain-containing protein 3, mitochondrial Alternative name(s): Mitochondrial transcription termination factor 2 Short name= mTERF2 Short name= mTERFL Mitochondrial transcription termination factor-like pr
Expression Region	36-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.