



# Recombinant Bovine DnaJ homolog subfamily B member 11 (DNAJB11)

<b>Product Code</b>	CSB-BP670243BO
<b>Abbreviation</b>	DNAJB11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q3ZBA6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GRDFYKIL GVPRSASIKD IKKAYRKLAL QLHPDRNPDD PRAQEKFQDL GAAEVLSDS EKRKQYDTYG EEGLKDGHS SHGDIFSHFF GDFGFMFGGT PRQQDRNIPR GSDIIVDLEV TLEEVYAGNF VEVVRNKPVA RQAPGKRKCN CRQEMRTTQL GPGRFQMTQE VVCDECPNVK LVNEERTLEV EIEPGVRDGM EYFIGEGEP HVDGEPGDLR FRIKVVKHSI FERRGDDLYT NVTISLVESL VGFDMDITHL DGHKVVHISR KITRPGAKLW KKGEGLPNFD NNNIKGSLII TFDVDFPKEQ LSEEAREGIK QLLKQGSVQK VYNGLQGY
<b>Source</b>	Baculovirus
<b>Target Names</b>	DNAJB11
<b>Protein Names</b>	Recommended name: DnaJ homolog subfamily B member 11 Alternative name(s): ER-associated DNAJ ER-associated Hsp40 co-chaperone ER-associated dnaJ protein 3 Short name= ERdj3 Short name= ERj3p
<b>Expression Region</b>	23-358
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.