



Recombinant Human Uncharacterized protein C1orf177 (C1orf177)

Product Code	CSB-EP670268HU-B
Abbreviation	C1orf177
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3ZCV2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRESQDAAGA HGWNRVGSTA TKWFTGAPFG VQSHRFDISA VYPNWKKFST FTEAPYSTRY STQVSHIGPG TYSSKETCFS KKKLMKEVDT GWAQAQEATR LTQLPHFQYQ AIMKEKRLKE QKLGPGSYNL KDFLEQLREK PCSTRGLLSS GEVRFRLTG NYYPGPGNYG EKGPNPYTKLE ENAWNRSRSHSE GLMCRMSNKP HRPYQGSGL GPGTYFFKSD LETYVARSVG TRGPYDTFSG DRSKPLPYGH YSMQKKKPRE LMNFKSFVEE LNSHHNKKKHG VFSKLPRNPK TPTERIYWAN LSQCPRTLAT SGPSFWLPQE KKCKPVNQPP FLLTSKGS KACQMIMGSW NPVGVGRYLN TWLMETKDRR QRYRSLFLSG SKRYLSDLAR DMLMQRITP FTKGKCPPTV DYNSDPTP
Source	E.coli
Target Names	LEXM
Protein Names	Recommended name: Uncharacterized protein C1orf177
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.