



Recombinant Human Paired box protein Pax-6 (PAX6)

Product Code	CSB-BP017492HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P26367
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQNSHSGVNV LGGVFNVRP LPDSTRQKIV ELAHSGARPC DISRILQVSN GCVSKILGRY YETGSIRPRA IGGSKPRVAT PEVVSIAQY KRECPSIFAW EIRDRLLESEG VCTNDNIPSV SSINRVLRLN ASEKQQMGAD GMYDKLRMLN GQTGSWGTRP GWYPGTSVPG QPTQDGCQQQ EGGGENTNSI SSGEDSDEA QMRLQLKRKL QRNRTSFTQE QIEALEKEFE RTHYPDVFA RERLAAKIDLPEARIQVWFSN RRAKWRREEK LRNQRQASN TPSHIPSS FSTSVYQPIP QPTTPVSSFT SGSMLGRTDT ALTNTYSALP PMPSTMAN LPMQPPVPSQ TSSYSCMLPT SPSVNGRSYD TYTPPHMQTH MNSQPMGTSG TTSTGLISPG VSPVQVPGS EPDMSQYWPR LQ</p>
Source	Baculovirus
Target Names	PAX6
Protein Names	Recommended name: Paired box protein Pax-6 Alternative name(s): Aniridia type II protein Oculorhombin
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes paired box gene 6, one of many human homologs of the Drosophila melanogaster gene prd. In addition to the hallmark feature of this gene family, a conserved paired box domain, the encoded protein also contains a homeo box domain. Both domains are known to bind DNA, and function as regulators of gene transcription. This gene is expressed in the developing nervous system, and in developing eyes. Mutations in this gene are known to cause ocular disorders such as aniridia and Peter s anomaly. Alternatively spliced transcript variants encoding either the same or different isoform have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.