



Recombinant Escherichia coli Uncharacterized protein yqeB (yqeB)

Product Code	CSB-BP670501ENV
Abbreviation	yqeB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q46808
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MNIFTEAAKL EEQNCPFAMA QIVDSRGSTP RHSAQMLVRA DGSIVGTIGG GMVERKVIEE SLQALQERKP RLFHGRMARN GADAVGSDCG GAMSVFISVH GMRPRLVLIG AGHVNRAIAQ SAALLGFDIA VADIYRESLN PELFPPSTTL LHAESFGAAV EALDIRPDNF VLIATNNQDR EALDKLIEQP IAWLGLLASR RKVQLFLRQL REKGVAAEEHI ARLHAPVGYN IGAETPQEIA ISVLAIEILQV KNNAPGGLMM KPSHPSGHQL VVIRGAGDIA SGVALRLYHA GFKVIMLEVE KPTVIRCTVA FAQAVFDGEM TVEGVATARLA TSSAEAMKLT ERGFIPVMVD PACSLLELK PLCVVDAILA KQNLGTRADM APVTIALGPG FTAGKDCHAV IETNRGHWLG QVIYSGCAQE NTGVPGNIMG HTRRVIRAP AAGIMRSNVK LGDLVKEGDV IAWIGEHEIK APLTGMVRGL LNDGLAVVGG FKIGDIDPRG ETADFTSVSD KARAIGGGVL EALMMLMHQG VKATKEVLEV A
Source	Baculovirus
Target Names	yqeB
Protein Names	Recommended name: Uncharacterized protein yqeB
Expression Region	1-541
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.