



# Recombinant Escherichia coli Protease elaD (elaD)

<b>Product Code</b>	CSB-MP670613ENV
<b>Abbreviation</b>	elaD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q47013
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMVTVVSNYC QLSQTQLSQT FAEKFTVTEE LLQSLKKTAL SGDEESI ELL HNIALGYDKF GKEAEDILYH IVRTPTNETL SIIRLIKNAC LKLYNLAHIA TNSPLKSHDS DLLLFFKKLFS PSKLMTIIGD EIPLISEKQS LSKVLLNDEN NELSDGTNFW DKNRQLTTDE IACYLQKIAA NAKNTQVNYP TGLYVPYSTR THLEDALNEN IKSDPSWPNE VQLFPINTGG HWILVSLQKI VNKKNKLQI KCVIFNSLRA LGYDKENSLK RVINSFNSEL MGEMSNNNIK VHLNEPEIIF LHADLQQYLS QSCGAFVCMA AQEVIEQRES NSDSAPYTLL KNHADRFKKY SAEEQYEIDF QHRLANRNCY LDKYGDANIN HYYRNLEIKH SQPKNRASGK RVS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	elaD
<b>Protein Names</b>	Recommended name: Protease elaD EC= 3.4.22.- Alternative name(s): Deubiquitinase Deubiquitinating enzyme Short name= DUB Deubiquitinating protease
<b>Expression Region</b>	1-403
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.