



Recombinant Human Tripartite motif-containing protein 60 (TRIM60)

Product Code	CSB-MP670964HU
Abbreviation	TRIM60
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q495X7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEFVTALVNL QEESSCPICL EYLKDPVTIN CGHNFCRSCL SVSWKDLDDT FPCPVCRCFC PYKSFRNPQ LRNLTEIAKQ LQIRRSKRKR QKENAMCEKH NQFLTLFCVK DLEILCTQCS FSTKHQKHYI CPIKKAASYH REILEGSLEP LRNNIERVEK VILQGSKSV ELKKKVEYKR EEINSEFEQI RLFLQNEQEM ILRQIQDEEM NILAKLNENL VELSDYVSTL KHLREVEGK SVQSNLELLT QAKSMHHKYQ NLKCELFSS RLTKYGFSLP PQYSGLDRII KPFQVDVILD LNTAHPQLLV SEDRKAVRYE RKKRNICYDP RRFYVCPAVL GSQRFSSGRH YWEVEVGNKP KWILGVCQDC LLRNWQDQPS VLGGFWAIGR YMKSGYVASG PKTTQLLPVV KPSKIGIFLD YELGDLSEFYN MNDRSILYTF NDCFTEAVWP YFYTGTDSEP LKICSVSDSE R
Source	Mammalian cell
Target Names	TRIM60
Protein Names	Recommended name: Tripartite motif-containing protein 60 Alternative name(s): RING finger protein 129 RING finger protein 33
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein contains a RING finger domain, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.