



Recombinant Flap endonuclease 1 (FEN1)

Product Code	CSB-BP671198LOG
Abbreviation	FEN1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q4FYU7
Product Type	Recombinant Protein
Immunogen Species	Leishmania major
Purity	>85% (SDS-PAGE)
Sequence	MGILGLSKLL YDKSPNAIRE QELKNFFGRR IAVDASMSIY QFIAMKGFQ DGQGLELTNE KGDVTSHLNG LFARTLRMID EGIKPIYVFD GKPPKLKADE LEMRRQKAAE AERAFEKAKD AGDDEMMEKM SKRTVRVSRD QIDESKLLR LMGIPVIQAP SEAEAQCAEL VKKGKAWAVG TEDMDALTFG STVMLRHLNI SDAKKRPIVE IHLDEVLQTT GLSMDQFVDL CILLGCDYVP KVPGIGPQKA WEGIQRYGSI ESFLESDDTT KHPVPADFYY KEARAFFQNP EVTRAEEIDI QFSEPDEVGL IQFLVKEKLF NPDRVNKGIA RLRAAFTRKT QGR LDSFFTI TKVPQQTAAA RAPLVGTKRP RDGKYVHVSG TLRKATSGHK KAVKK
Source	Baculovirus
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>