



Recombinant *Oryza sativa* subsp. *japonica* Expansin-B1 (EXPB1a)

Product Code	CSB-EP671430OFG
Abbreviation	EXPB1a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q40638
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	GPPKVPPG PNITTSYGDK WLEAKATWYG APKGAGPKDN GGACGYKDVD KAPFLGMNSC GNDPIFKDGK GCGSCFEIKC SKPEACSDKP ALIHVTDMND EPIAAYHFDL SGLAFGAMAK DGKDEELRKA GIIDTQFRRV KCKYPADTKI TFHIEKASNP NYLALLVKYV AGDGDVVEVE IKEKGSEEWK ALKESWGAIW RIDTPKPLKG PFSVRVTTEG GEKIIAEDAI PDGWKADSVY KSNVQAK
Source	E.coli
Target Names	EXPB1a
Protein Names	Recommended name: Expansin-B1 Alternative name(s): Beta-expansin-1 Major pollen allergen Ory s 1 Ory s I OsEXPB1 OsaEXPb1.2 OsaEXPb1.3 Allergen= Ory s 1
Expression Region	23-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.